

Mastereo SRS SYSTEM



Model No. : R00SRS31-00

Panacea's Im.PAc.T Mastereo SRS system is designed with patient's comfort in mind and fulfils the requirement of achieving highest accuracy with daily reproducibility during Stereotactic Radiosurgery.



The carbon fiber base plate of the Mastereo SRS system is extremely light weight and can be indexed directly on any couch top. The lower end of the base plate is tapered to provide relaxed positioning to the patient's shoulder.

The unique head rest of the Mastereo SRS System is made of "transparent soft PVC" which conforms to the shape of the patient's head. This allows customized positioning for all the patients.



The Mastereo SRS mask provides high definition to the shape of the head. With instant hardening and lesser shrinkage, patient positioning time is reduced. The Mastereo SRS masks comes with different perforation size and pattern to blend with the patient's anatomy and exerts less pressure on areas which does not receive radiation, thereby ensuring right immobilization and comfort for patient.