

## “ For Validation, Correction & Registration”


Image Registration Phantom from Panacea is a distinctively designed platform for a three-dimensional registration of diagnostic imaging data with the planning computed tomography which offers a substantial improvement in tumor target identification for many radiation therapy patients. The purpose of this phantom serves the 6D couch shift corrections and validations.

Insert Material	Tissue Equivalent & Bone Equivalent
Insert Shapes	Irregular, Circular & Oval
Phantom Layers	5 layers for motion & 1 fixed layer
Motion Range	-5° to +5° in 1° resolution

Model No  
R00IRP07-00



Rigid & Meticulous construction



Accurate



Quick Setup



Fine tune adjustments

